

TLi - PRO POWER C BOOSTING MASK

This LIGHTWAVE Topical Light Infusion® (TLi) Only vitamin C treatment is a two-step process and is 12 minutes in duration. Adding **Lightwave Red Light Therapy*** will increase the treatment time. A highly concentrated Power Pure™ freeze-dried Vitamin C sphere is blended with scientifically developed Peptides, Sepitonic™ M3 Active Mineral Complex and Glutathione for a dose of triple-action, lightening and brightening and age fighting super-stars. This innovative development in skin care allows for each topical light infusion treatment to remain fresh at its peak potency until activated for use. This 12-minute treatment calms the skin while dosing it with a potent amount of lightening and brightening ingredients for a more beautiful, radiant, and youthful appearance.

Note: To activate the Topical Light Infusion® setting on the system, you must first purchase Topical Light Infusion®. Only Topical Light Infusion can be used with this system. Additionally, take all necessary before pictures of target area prior to starting session number one. If the photos are taken at the conclusion of the treatment, the skin can temporarily appear flushed and red in color.

1. Prepare the Skin: It is important to properly prep the skin for Topical Light Infusion®. This is done by thoroughly cleansing the skin with the abi's Corrective Wash. This will ensure all excess oil, make-up, dirt, and debris has been removed. For significantly aged and/or sun damaged skin, additional exfoliation is recommended prior to starting the Topical Light Infusion® treatment. Microdermabrasion or abi's various scrubs and chemical peels are all highly effective in properly preparing the skin.

**Lightwave LED Light therapy option: Including an LED light therapy session with your Topical Light Infusion® treatment will provide the best possible results. This can be done by completing a standard LIGHTWAVE treatment before completing the Topical Light Infusion® program. If you are limited on time, you can pre-treat the skin by turning on the "LW" button at the bottom of your screen prior to starting the TLi program. This will add an IR ONLY LED light session to the beginning of your treatment and increase your treatment time.*

2. Topical Application: Step 1 of 2, remove cap from the Pro Power C Sphere Jar. Crush the Pro Power C ball in the glass jar using the spatula found on the lid. Remove the cap from the Step 1 Liquid Activator syringe and disperse the liquid into the jar containing the crushed sphere. Replace the lid on the jar and shake well. Use the syringe to draw the liquid back out of the jar. Apply the entire liquid to the target treatment area. You have been provided with enough serum to cover the face and neck area. Allow serum to penetrate the skin for 30 – 60 seconds before moving on to step 2.

Next, Step 2 of 2, remove the Pro Power C Gel Mask from foil packet. Gently unfold the topically infused gel mask.

Remove the paper backing and place the mask on the face with the eyes, nose and mouth openings positioned in the appropriate places. To maximize the mask, smooth the gel so that it lays flat on the skin and covers as much area as possible but still allows the client to feel comfortable. Apply any excess liquid from the packet to the décolleté area.

Precautions: Once both serum and mask has been applied, shield the patient's eyes by fitting the protective goggles in place. The use of goggles is required to minimize brightness and to avoid any incidental eye exposure.

3. Treatment Instructions Place the LED panel directly over the face and neck area. To maximize the results of the treatment, make sure all six panels are as close to the target area as comfortably possible. Activate the TLi treatment. The LED dome panel will remain lit with red light for 12 minutes. The light is pulsing even though it appears to be constant. Therefore, individuals with epilepsy should be monitored closely.

4. Post Care: Finally, at the end of the Topical Light Infusion treatment, remove the LED panel. Next, apply a dime size amount of the ABI Intense Antioxidant to help nourish the skin and dispose of free radicals. Allow the product to penetrate. Then apply a dime size amount of abi's Hydrating Recovery Gel and abi's Moisturizing Antioxidant to help protect the skin's natural moisture factor. Finally, apply sunscreen to protect the skin from any future UV damage. It is also important to always have the patient increase their water intake by at least 8oz following a treatment.

NUMBER OF SUGGESTED TREATMENTS: This mask is used as a stand-alone 1-4 session treatment or as a Topical Light Infusion series. When using it as a stand-alone series, it is recommended not to do this mask more than twice a week. If it is used in a series, start the series with this mask and follow with 3 Wrinkle Repair TLi treatments and then insert one to two more boosting gel masks throughout the series. For maintenance, it is recommended to perform a TLi treatment every 4-6 weeks or as needed and apply the Pro Platinum C at home daily. The Combination of LIGHTWAVE™ and Topical Light Infusion® are highly recommended for the best possible outcome.